




ONLINE STORE
Purchase@SYntonac

Product Catalog


Experts make things simpler



SYntonac International LLC

office@syntonac.com

www.syntonac.com



SYntonac International LLC

was founded on Syntonic Optometry principles as described by H.R. Spitler, and S.W. Collier.

*"Syntony" Being adapted by a certain wavelength.
Quote by H.R. Spitler.*

The Company directors are: Drs. Vasana Lerdvoratavee, Chula Lerdvoratavee and Ronald Weingart.

SYntonac International is headquartered in Oswego, IL - USA.

We produce a quality line of products that allows the Behavioral Syntonic Optometrist to use a variety of

specific filters and frames to improve the therapeutic outcomes of their patients.

We also provide education for treatment protocols for each mode of therapy.

We will be present at several professional events, where you can discover our New SYntonac Line of products.

Visit the Seminar- and Event section on our website to find out where you can meet us.

Want to know more about our products?

Please contact us, we're always happy to answer any questions.

2075 Wiesbrook Road,
Oswego, IL 60543
USA

SYntonac Office - USA



office@syntonac.com

www.syntonac.com

4

Focal - Local

SYntonac Penlight Stimulation for the Optometrist needing precision in their Treatment

6

Filter Goggles

Ideal for in-office and at home treatment regimens. These are based on Syntonic principles where the first filter in a combination of filters is closest to the eye.

8

Syntonizer

Ideal for in-office and at home treatment regimens. We have 2 Syntonizers. The "ELAN" and "EASY" Units. These can be used via 2 different Modalities.

10

Field Charter

The SYntonac Field Charter serve many functions. It is able to discover and demonstrate functional vision deficits. And much more...

11

Monocle Paddle

Use for: Retinoscopy Filter Testing and/or monocular procedures. The Monocle Paddles are included with the "ELAN" Syntonizer Unit..

12

Binocle Paddle

Use for: Testing, Binocular procedures and/or the Asymmetric Equilibrium Procedure (created by Dr. Chula Lerdvoratavee).

13

Clear Slant Board

The use of the SYntonac Slant Boards is based on the educational work by Darrell Boyd Harmon, PhD.

13

Clear Waddle Board

The Waddle Board is a Multi-functional Balance Board.

14

FSO Forum

The FSO Forum is supported by a team collaboration software that brings people and the resources they need to one place, making it easier to share and learn about FSO.

15

Accessories

We have items that can improve or deepen the therapy. And we also have items that improve the organization of your office.

We offer products to treat the patient through three main modes.

Because of our Company policy, no product pricing will be displayed in this brochure.

About

Focal - Local Syntonic Optometry

- ✓ Is a method applying specific light wavelengths to extrinsic muscle points, as other vision points on the body to bring about an energetic balance to improve visual function.
- ✓ It has been tested on a variety of vision problems, with Syntonic Optometric testing before and after to determine the immediate response.



We offer two Focal - Local Syntonic Penlights:

- *In-office case*
- *Home unit*



- ✓ We offer a wide array of filter combinations.

- ✓ The SYntonac Local-Focal line is developed by Yvonne Frei and Stefan Collier.

The Focal - Local Penlights are also ideal to use as:



- ✓ Syntonization Filter Penlight **with** or without the "EASY" Syntonicizer unit.
- ✓ Syntonization during VT Procedures such as:
 - Penlight Pursuits
 - Penlight Tromboning
 - Penlight Push-up
 - Penlight Push-away
 - Penlight Saccades
 - Penlight Fusion
 - Clear Maddox with Penlight Procedure

SYntonac Office Set

The SYntonac-Office set is the ideal Kit for the optometrist to use in the Office.

Ideal



Office

Set includes

- ✓ 22 Syntonac Penlights, 1 white Maddox paddle, 2 fixation white VT-penlight, 12 Foxy/Grey frames with the 12 basic Syntonic filter combinations, 1 case.
- ✓ The 21 SYntonac Penlight filter combinations are:
 - Alpha, Alpha-Delta, Alpha-Theta, Alpha-Omega, Alpha-Lambda, Depressant, Stimulant, Mu, Mu-Delta, Mu-Theta, Mu-Upsilon, Mu-Pi, Omega, Upsilon, Pi, Upsilon-Omega, Upsilon-Omega-D, Pi-Omega, Theta, Neurasthenic, Omega-N.



SYntonac Home Unit

Ideal Home Unit for Local-Focal Syntonic Optometry.

Ideal



Home Office

Syntonization Filter Penlight to use with or without the "EASY" Syntonizer Unit.

Set includes

- ✓ Comes in a box with 1 penlight and 1 battery.
- ✓ You can choose your Filter Combination.
- ✓ We offer a wide array of filter combinations.



SYntonac Clear Unit

Ideal Home Unit for Local-Focal Syntonic Optometry.

And to use with or without the "EASY" Syntonizer Unit.

Ideal



Home Office

Syntonization Clear Penlight to use with or without the "EASY" Syntonizer Unit together with SYntonac Filter Goggles.

Set includes

- ✓ Comes in a box with 1 clear penlight and 1 battery.



Silver Fixation Penlight

The SYntonac Silver Penlight will provide a more integrative result.

It is ideal for Home and Office Syntonic Therapy.

Ideal



Home Office

This silver non-LED penlight is an ideal light source and fixation target to be used together with the SYntonac Goggles.

Set includes

- ✓ 1 Penlight
- ✓ Battery included



About

SYntonac Filtered Goggles

- ✓ We offer a wide variety of filter combinations to meet the specific needs of each individual patient.
- ✓ We also produce Nascentization goggles for Local and/or Non Local stimulation.
- ✓ We offer a truly precise product.



The Filter Goggles are available in:

- ✓ Basic Sets
- ✓ Per Piece
- ✓ Customized Combo
- ✓ Six Pack
- ✓ Discover All Possibilities and Filter Variations in the Online Store

- Office
- Home



Frames for Adults:

- ✓ We offer two frame designs:
 - Foxy
 - Boxy
- ✓ **Foxy** is our stylish general filter goggle. It allows a wider peripheral stimulation during therapy. It is a wonderful design for most treatment protocols.
- ✓ **Boxy** is our specialty goggle specifically designed with a therapeutic component. It is designed to control the amount of light entering the eyes, helping to rebuild the structure in the injured visual system.

Frames for Children:

- ✓ We offer two frame designs:
 - Maxy
 - Mixy
- ✓ **Maxy** is our stylish general filter goggle. It allows a wider peripheral stimulation during therapy. It is a wonderful design for most treatment protocols.
- ✓ **Mixy** is our specialty goggle specifically designed with a therapeutic component. It is designed to control the amount of light entering the eyes, helping to rebuild the structure in the injured visual system.

Frame Colors

We offer the Foxy-Boxy and Maxy-Mixy frame designs in three different color frames:

- ✓ White
- ✓ Grey
- ✓ Black

Why three different color frames?

Create your Custom Multiple Filter Goggles

You have the possibility to create your own Filter combination up to 4 filters:

- ✓ You select the Basic Filters and the next layers.
- ✓ The Basic Filter is the first filter closest to the eye. You can add up to 3 filters to combine with the first basic filter.
- ✓ The Basic Filters are: Alpha, Delta, Theta, Mu, Omega, Upsilon, Pi, Lambda

Many filters combinations possible. - Select the goggle frame and color and the filter layers.

Specialized Filter Goggles / Special Combo

Since our society has become physically, emotionally, and visually more complex and demanding, our patients/clients need to be approached in a more precise and functional manner.

Functional Syntonic Optometry (FSO) adds another, deeper dimension to current treatment modalities.

It fulfills treatment demands with more advanced measurements, different and a greater variety of training techniques adapted to each patient's personal needs.

Therefore we have created Special Filter combinations adapted to the patient's personal needs.

- ✓ The goggles will be applied (in most cases) with the FSO Trans-lid technique.
- ✓ Many filter combinations possible.

Stefan Collier has discovered that White, Grey and Black frames have a different influence on Magno cellular processing during Syntonic Optometric therapy.

- ✓ The white frame gives a more stimulating impulse to the Magno cellular processing.
- ✓ The grey frame gives a more balanced - neutral impulse to the Magno cellular processing.
- ✓ The black frame gives a more relaxing impulse to the Magno cellular processing.

FILTER GOGGLES

We Offer

a wide variety of filter combinations to meet the specific needs of each individual patient.

- Carton goggles with Multiple Filter Layers
- Specialized Filter Goggles
- Nascentization Goggles
- Vision Therapy Goggles

Truly precise products

Nascentization Filter Goggles

We have Nascentization goggles for Local and/or Non Local stimulation.

Goggles for Vision Therapy

We have special VT Kits available :

- ✓ Magno cellular VT Kit
- ✓ Gross Motor Stretching Kit
- ✓ Eso Myopia Kit
- ✓ EXO Presbyopia Kit

These VT Kits are available in Adult or Children Frames.



About

SYntonac Syntonizer

We have 2 Syntonizer Units available:

- ✓ "ELAN" Unit
- ✓ "EASY" Unit

- Office
- Home



The "ELAN" Unit includes

- ✓ An additional spare Light Unit
- ✓ A set of 21 Monocle Paddles

The Monocle Paddles are:

- ➔ Alpha - Omega, Alpha - Delta, Alpha - Theta, Mu - Delta, Mu - Theta, Mu - Upsilon, Omega - N, Pi - Omega, Upsilon - Omega, Upsilon - Omega - D, Alpha, Delta, Theta, Mu, Omega, Upsilon, Pi, Lambda, D, N, S

Syntonizer "ELAN"

This instrument is designed for in-office treatment regimens.

The "ELAN" Syntonizer Unit is:


- ✓ Light Weight and Slightly Larger than Easy Unit.
- ✓ Very Durable, Dimensionally Stable and non-toxic synthetic material.
- ✓ A Light Source at the far end of the instrument in front of four Monocle Paddle Slots, allowing for a multitude of Filter combinations.
- ✓ The Light Source has a USB plug, allowing Worldwide access.
- ✓ The "ELAN" Unit includes an additional spare Light Unit to ensure no lost Therapy Time.
- ✓ Versatile: Making it possible to Syntonize in two different Modalities
 1. The instrument can be operated with SYntonac Filter Goggles.
 2. Or any combination of Filters inserted into the Monocle Paddle Slots.
- ✓ Designed with integrated Fixation Lines, allowing the patient to engage in Ocular Motor Vision Training while Syntonizing.
- ✓ Adjustable in different height positions for patient comfort.

Syntonizer "EASY"

This instrument Ideal for in-office and at home treatment regimens. The "EASY" Syntonizer Unit is:

The "EASY" Syntonizer Unit is:

- ✓ Very light.
- ✓ Very durable, dimensionally stable and non-toxic synthetic material.
- ✓ Versatile: Making it possible to Syntonize in two different Modalities
 - ▶ The instrument can operated with the clear SYntonac Penlight together with the SYntonac Filter Goggles.
 - ▶ Or choose the Focal-Local Syntonic Penlight combination your patient needs.
- ✓ Designed with integrated fixation lines, allowing the patient Vision Training while Syntonizing.
- ✓ Adjustable in different height positions for patient comfort.



- ✓ Discover All Possibilities Syntonizer Variations in the Online Store
- ✓ We created a special "EASY UNIT" Carry-on Bag. This bag gives the ideal protection for the "EASY UNIT", Filter Goggles and Filter Penlights. Ideal for Home Therapy use. The bag is sold separately.

How to use:

Modality 1

With the SYntonac Clear Local-Focal Penlight and SYntonac Filtered Goggles:

1. Remove the insert on the back of the "EASY" Syntonizer.
2. Turn on the penlight.
3. Place the **Clear Local-Focal penlight** in the therefore designed hole.
4. Replace the insert into the back of the "EASY" Syntonizer
5. Put on the selected SYntonac Filter Goggles.
6. Carefully move your face into the opening of "EASY" Syntonizer
7. You're ready to start the Training and syntonization

How to use:

Modality 2

With the SYntonac Local-Focal Penlight with the integrated Filter Frequency:

1. Remove the insert on the back of the "EASY" Syntonizer.
2. Turn on the penlight.
3. Place the SYntonac Filtered Local-Focal penlight in the therefore designed hole.
4. Replace the insert into the back of the "EASY" Syntonizer
5. Carefully move your face into the opening of "EASY" Syntonizer
6. You're ready to start the Training and syntonization

About

SYntonac Color Field Charter

Color Field Charters are the newest addition to our product line.

- ✓ They are fabricated from durable, light-weight acrylic polymers in two colors - Clear and Black.
- ✓ Each unit is illuminated with a light source that has three settings.
- ✓ A moveable camera is able to view each eye during testing.
- ✓ Two complete sets of testing probes are included. Each set has five target sizes in each color - Red, Blue, Green and White.



• Office




The "SCFC" Unit

Field Charters serve many functions.

- ✓ to discover and demonstrate functional vision deficits.
- ✓ to determine if the Syntonac Therapy is working.
- ✓ assess the performance of the Filter frequency and strength

Available in Clear and Black:

- ➔ The **Clear Field Charter** is not only beautiful and top design, but it gives the patient the sensation of looking in open space. Which results in a more natural measurement



This affordable Color Field Charter is a welcome addition to any Syntonac practitioner..

SYntonac Monocle Paddle

The Monocle Paddle in Functional Syntonic Optometry is used for:

- ✓ Retinoscopy Filter Testing
- ✓ Monocular procedures
- ✓ Syntonizer "ELAN" Unit



The Monocle Filter Paddles are available in:

- ✓ Different Color Frames:
 - ▶ Black,
 - ▶ White,
 - ▶ Grey
 - ✓ Basic Sets
 - ✓ Customized Combo
 - ✓ Per Piece
 - ✓ Six Pack
- ✓ Discover All Possibilities and Filter Variations in the Online Store

- Office
- Home



How to use the Monocle Paddle

Directions for using SYntonac Monocle Paddle:

1. Filter Combinations - There will be a different Filter Combination on each side of the paddle.
2. The Filter Combination facing the patient is the Filter that the patient sees.
3. That is, if the Filter Combination that faces the patient reads Omega-N, Omega is the filter closest to the the patient and N is behind it and further away from the patient's eye.
4. The optometrist sees the opposite designation of the paddle.



Retinoscopy Filter Testing



Monocular procedures



Syntonizer "ELAN" Unit

The Binocle Filter Paddles are available in:

Available as

Basic Sets - Per piece - Six Pack - Customized Combo

Different Color Frames:

Black - White - Grey

Discover

All Possibilities and Filter Variations in the Online Store

- Office
- Home



SYntonac Binocle Paddle

Use for:

- ✓ Testing
- ✓ Binocular procedures
- ✓ Asymmetric Equilibrium Procedure (created by Dr. Chula Lerdvoratavee)



SYntonac Waddle Board

The Clear Waddle Board is a Multi-functional Balance Board.

The Benefits of a Clear Waddle Boards are the following:

- ✓ Stimulates and Balance Peripheral Awareness
- ✓ Supports Multi-tasking
- ✓ Enhances Eye/Hand and Eye/Foot Coordination
- ✓ Promotes Figure-Ground Relationships
- ✓ Fortifies Body Awareness
- ✓ Reinforce Balance

Our Clear Waddle Boards are available in two variants:

- ✓ Short: 18" x 6" x 4"
- ✓ Long: 24" x 6" x 4"

Both variants are Durable, Light and Portable.

- Office
- Home



SYntonac Slant Desk

The application of the Clear Slant Board is every important for your Vision.

The use of the SYntonac Slant Boards is based on the educational work by Darrell Boyd Harmon, PhD.

The Benefits of a Clear Slant Boards are the following:

- ✓ Precisely designed at ideal 20 degree slant
- ✓ Reduces Reading Stress by 40%
- ✓ Promotes Figure-Ground Relationships
- ✓ Stimulates and Balances Peripheral Awareness
- ✓ Complements Vision Training Procedures
- ✓ Enhances Near and Global Reading

Our Clear Slant Boards are available in two variants:

- ✓ Small: 14" x 11" x 6"
- ✓ Large: 22" x 14" x 6"

Both variants are Durable, Light and Portable.

Offering a precise product leads to a faster and more successful therapy.



The FSO Forum

is supported by a team collaboration software that brings people and the resources they need to one place, making it easier to share and learn about Functional Syntonic Optometry (FSO).

Get together to get more experienced.

There're 16 different Topic groups to give each subpart of Functional Syntonic Optometry an exclusive space where sharing can happen without barriers.

Communicate faster and collaborate better.


Keep in touch with your colleagues, no matter where you are, using Feeds. Exchange ideas, share

cases, ask for help, learn through reading the FSO Support Articles, receive updates from FSO Forum members, and always stay up to speed on what's happening in FSO.


Keep FSO files where you can find them.

Bring FSO files to one place, so all of you can access them when needed. Create documents with your colleagues using the online editors and turn writing into a collaborative process.

- ✓ You, as a Functional Syntonic Optometrist, but most of all, your patients will have great benefit from it.



✓ Offering a precise product leads to a faster and more successful therapy.



A Accessories

We have items that can improve or deepen the therapy. And we also have items that improve the organization of your office.

I Carry-on Bag "EASY UNIT"

We created a special Carry-on Bag for our "EASY UNIT" SYntonizer.

This bag gives the ideal protection, store and transport for the "EASY UNIT", Filter Goggles and Filter Penlights.

✓ Ideal for Home Therapy use.

2 Penlight Holder

Bringing your favorite filter penlights together in one place.

✓ Fits 12 penlights

✓ The holder is fabricated from durable, light-weight acrylic polymers.

3 Clear Maddox

This Clear Maddox is ideal to use together with the SYntonac Goggles and or SYntonac Penlights to obtain a more integrative result. Ideal for Home and Office therapy.


This Procedure with the clear maddox is for integration of Phorias.

✓ Such as: Eso, Exo, Hyper, Hypo Phoria. You can use it for near and far.



Sales Conditions

Experts make things simpler



Syntonac International LLC

office@syntonac.com

www.syntonac.com

Welcome to the SYntonac Online Shop

Our Company policy is that we only sell to Optometrists - Optometry Offices.

- ✓ Therefore if you want to login in our SYntonac Online shop you need a Password.
- ✓ Fill out the Online Login Password Request Form.
- ✓ After Profession Identification check you will receive an Email with your Password.



Buy Direct

Why buy direct?

- ✓ We save you time and capital, while still delivering the finest Syntonic Optometry equipment.

30-day money back guarantee

At SYntonac International, we want you to be 100% satisfied with our SYntonac Line.


- ✓ In the unlikely event that you are not completely satisfied with your purchase, you may return unopened items for a full refund within 30-days of delivery.

Shipping

Worldwide shipping

- ✓ Affordable, worldwide shipping available. Receive a shipping quote on our checkout page!
- ✓ Shipping is based on parts availability.
- ✓ Please allow 3 business days for processing for all orders.

Feel free to contact us to discuss shipping charges, questions, or concerns.



SYntonac International Light System

SYntonac

Catalog